

3M™ ESPE™ Strip Crown Forms – Pedo | Size Overview

Upper central incisor (Patient's left)	Upper central incisor (Patient's right)	Upper lateral incisor (Patient's left)	Upper lateral incisor (Patient's right)
 611 UL1 6.0 mm 914001	 511 UR1 6.0 mm 914011	 621 UL1 4.3 mm 914021	 521 UR1 4.3 mm 914031
 612 UL2 6.7 mm 914002	 512 UR2 6.7 mm 914012	 622 UL2 5.1 mm 914022	 522 UR2 5.1 mm 914032
 613 UL3 7.4 mm 914003	 513 UR3 7.4 mm 914013	 623 UL3 5.9 mm 914023	 523 UR3 5.9 mm 914033
 614 UL4 8.1 mm 914004	 514 UR4 8.1 mm 914014	 624 UL4 6.7 mm 914024	 524 UR4 6.7 mm 914034

